

# All in one time tagging, coincidence correlation and delay/pulse generation

Observe detection rates for up to five single-photon detectors, generate up to four coincidence histograms between any pair of connected detectors or trigger signals, and record every time tag for an exceptionally high number of detected photons, all with picosecond precision.

Beyond measurement, access four outputs with customizable and conditional pulses, to control and react to your experiments and applications in real time.

Get the most out of your photonic experiments and applications with the ID1000 Time Controller Series today.

## APPLICATIONS

- QKD and quantum communication
- Quantum optics and computing
- Single-photon source characterisation
- ► Fluorescence lifetime imaging
- ► Failure analysis of integrated circuits
- ► VIS, NIR and MIR spectroscopy

## **KEY FEATURES**

#### Five interchangeable input and reference channels

- High-resolution mode: up to 300 MHz across all inputs, with 1 ps resolution
- High speed mode: up to 100 MHz per input, with 100 ps resolution
- Measure four coincidence histograms in parallel
- **Record the precise arrival time of every detection event**, with or without coincidence filtering
- Tailored pulses (NIM or TTL) from four output channels, with customizable patterns and widths
- Advanced integrated logic: filter events for two-, three- and four-fold coincidences in real time
- **Quick and easy lab integration:** suite of LabView and Python scripts

## NEW

- Measure more: multi-device synchronisation for over 64 input channels
- Enhanced precision: 1 ps resolution with < 4 ps rms jitter, and excellent DNL for low signal applications

## **ID1000 Time Controller Series**

#### Counting and histograms

- Count single channel input rates
- Count 2/3/4-channel correlation rates
- Record four histograms simultaneously
- Up to 16,384 time bins, each between 1 ps and 4 ms wide

#### **Time tagging**

- Record every detection event, up to four timestamp measurements simultaneously
- Realtime 2/3/4-coincidence filtering
  Stream timestamps (raw or filtered) up
- to 10 M events/s

#### Output/delay generation

- Up to 4 x 250 MHz output pulses
- Customizable periodic or conditional on demand pulse sequences
- Digital NIM or TTL pulses, ideal for experimental instrument control



#### The last word in precision

Putting it all together

• 1 ps digital resolution (HR)

• Four programmable outputs

• Process up to 100 MHz per input

• < 4 ps rms jitter (single-channel, HR)

BO

Thanks to the advanced discrimination electronics of the ID1000 Time Controller, the precision and throughput of device will be future-proof for your lab, all in a versatile and cost-effective solution.

• Five interchangeable input and reference channels

(Right) Timing jitter characterization of an example ID1000 Time Controller Series device in high resolution mode. Measured as the root-mean-squared (rms) width of a start-stop histogram between an Input channel and Reference channel pair, with fast synchronized pulses, divided by a factor of  $\sqrt{2}$ .

## How to: TCSPC and the Instrument Response Function

The Instrument Response Function (IRF) describes the temporal response of your detection system against the 'true' dynamics of the system under investigation. The IRF takes into account all the detectors, measurement electronics, source precision, and any other instruments contributing to a non-instantaneous response in your experiments. Ultimately, the IRF tells you the best achievable timing precision of your detection system.

In the example to the right—typical for timecorrelated single-photon counting (TCSPC) the overall system response,  $R_{Sys}(\tau)$ , is a linear combination of the component IRFs, where  $\otimes$ is the convolution operator.

The resulting uncertainty – the timing jitter  $\Delta \tau_{Sys}$ -is thus the quadrature sum of all component timing jitters. Use this insight to find the best ways to improve your experimental precision.



(Above)  $\tau$ : time delay between detection events, 'Src': photon source, 'Det': singlephoton detector pair, 'TC': time correlator (e.g. ID1000), 'Rest': everything else in the detection system contributing to a finite timing response.

## ID1000: made for you

Simpler than ever to use and integrate into your laboratory, **ID1000 Time Controller Series** devices come with a new-and-improved suite of LabView virtual instruments and Python scripts.

**ID1000 Time Controllers** are available in six combinations depending on your needs (see Specifications on the next page), each upgradeable at any time. (Right) Example instrument response function of an ID281 SNSPD, at 1550 nm, recorded with an ID1000 Time Controller. The FWHM timing jitter of 24.4 ps includes all instrument jitter contributions, such that the real FWHM of the SNSPD is less than 22.3 ps.





(Above) Control, operate and display measurements of your ID1000 Time Controller with the included GUI application.

(Left) Directly access and intuitively configure the ID1000 Time Controller's internal logic, with the included Configuration Editor application, for limitless customization.

Bo

										Property	Value	
										HIST1:BCOU HIST1:BWID	6001 5	
										MIST1:DATA	N Plot	
1000	19085		19008	100	150062		mains .		01793	HIGTL: EXAB:LINE	REC	
NON	- 10	101		101	• 82		-			WTOTI STOP		
NIN	+903	NON	100	MONE.		-	100	10.00	1 1	DTOTATE DOD	-	~
nini	++++	100	#N0	noni	#1N0					RIDITIRIN		
	-		mount		THEFTER					HIST1:REP:FILT	RISI	
0012	+10	600		NONE.	•1A	1 100	100 C		04792	HIST1:REF:LINK	STAR	
nini	• HK	NON	101	NONE	•••	/ 7475	-	10000				
NOM	+++++	NON NON	• 80	NONE	•m0	/ 145					n	
NONE	+00	104	•00	NO4	· 400					HIST1:STOP:LINK	78005	
	79008		THOOLS		THEOLO							
60%	•••	600		MON	*** /	1625	+10					
NONE.	144	+ 604		404	···· /•	19030	+0104 +		+			
	100			1 100	and a	460	+EN8					
					7							
	75084	1.00	150002		150009		-					
	100		1	22	100	7805	+101					
100	The second			and a	-	79056	•934 •		* · · · *			
NEW	+00	105	eno	/ MOHE	+IND	ALC	•000					
00.44	10000	100	•m	NOVE			6445		OUTL			
NONE	+80	104	· ###	NONE	+100	534	-					
NING	+++++++++++++++++++++++++++++++++++++++	NON	++++ *	NONE		1835	1	al	1 1			
NIN	400	101	•m /	NONE.	*in0	~						
	-		mont		THEFT							
96.42	+10.	504	100	ND16	•FA	100	8492					
nini	*#K	NON		NONE	****	1000						
NINE	+000	100		NONE	-80	100	+5100					
NONE	+00	/ 101	•20	NOVE	*00							
	79007		THOOLS		TREEDS							
00.45	+18	/ site	#18	NONE		1925	+10					
NDM	192	• / NON		MORE	• 900	7876	4070F +					
-676		V 101		ACH4		INC	+048					

## In good company

- Operate multiple ID1000 Time Controller devices in unison
- Up to >64 synchronized input channels





(Above) Example of four ID1000 Time Controllers in a multi-device configuration, giving access to 16 interchangeable input and reference channels.

## You'll need some photons

The ID1000 Time Controller Series is designed especially for IDQ's range of single-photon detection solutions, from our versatile and easy-to-use semiconductor SPADs, to the game-changing performance of our SNSPDs.



## SPECIFICATIONS

Functional specifications	High Speed	High Resolution <sup>(1)</sup>	MASTER	TCSPC		
# Inputs channels per device		<	<b>~</b>			
# Input channels in multi-device mode (2)	>	<ul> <li>✓</li> </ul>				
Min. time bin width (digital resolution)	100 ps	1 ps	✓	<ul> <li>Image: A second s</li></ul>		
Timing jitter (single-channel, rms) (3)	< 28 ps	< 4 ps	<b>~</b>	×		
Max. DNL (rms)	< 10 ps	< 0.2 ps	✓	<b>~</b>		
Max. input count rate (4)	5 x 100 MHz	300 MHz	<ul> <li>✓</li> </ul>	<b>~</b>		
Max. timestamp rate (5)	10 M e	✓	<			
Input voltage range	-2 V t	✓	<b>~</b>			
Input channel delay	1 ps to 4 ms	✓	<			
# Output channels		✓				
Output pulse format	NIM	<ul> <li>✓</li> </ul>				
Max. rate per output channel	250	✓				
Min. output pulse width	6 ns (TTL),	✓				
Input — Output latency	400 ns	1500 ns	<ul> <li>✓</li> </ul>			
Electrical & Environmental specifications						
Power supply requirements	100 V to 240 V, 50 Hz to 60 Hz, 1 A to 2.5 A					
Operational temperature	5°C to 35°C					
Max. environment humidity	80% up to 31°C, 50% up to 35°C					
Device dimensions (W x H x L)	380 mm x 63 mm x 255 mm					

(1) HR add on license required.

(2) Number of accessible input channels when 16 Time Controllers are synchronised

(3) Measured as the width of a two channel start stop histogram with fast input pulses, divided by a factor of  $\sqrt{2}$  to give the single channel jitter. Timing jitter for high speed mode is a worst case estimate with an assumed Gaussian instrument response profile.

(4) Count rates available for single channel counters, coincidence counters, and histogram processing, across all input channels of a single device. (5) Measure up to four sets of timestamps (detection event time tags) in parallel, in total up to the max. timestamp rate.

#### On request, the two core ID1000 versions can be remotely upgraded with these add-on licenses:

Available add and	Core ID1000 version			
	MASTER	TCSPC		
HR - Access to 'high resolution' (1 ps resolution, < 4 ps jitter) mode for the 5 input channels	٤	٤		
<b>PRCSG</b> - Access to the internal FPGA logic, to configure real time selection filters for up to four fold coincidences	✓	٤		
4 OUT - Access to the 4 output channels and the internal delay / pulse generator functionality	<	٤		

🛃 Upgrade available

already included



## www.idquantique.com | info@idquantique.com

USA Boston, MA

South Korea SungNam-si